



Specifications

Model	SANDRIDER
Overall Length	72"
Overall Width	38"
Overall Height	48"
Assembled weight	49 pounds
Front wheel	11.8"
Rear wheels	19.3"
Recommended Maximum Weight	350 pounds
Seat height	19"
Seat width	18"
Seat depth	16"
Inter armrest spacing (standard)	16"
Inter armrest spacing (comfort fit)	21"

Materials:

Frame: Powder coated* aluminum 1.7" diameter frame, 1.3" diameter push handle and armrests

Upholstery*: Sling quality vinyl coated polyester, UV stable, Microban® antimicrobial protection, GREENGUARD® Select Certified

Fasteners: 316 stainless steel bolts and nuts with nylon washers

Wheels: Wheeliez™ low pressure polyurethane tires, polypropylene hub, 1" bearing

(*Custom upholstery fabrics and frame coating materials may vary)

patent pending